


CAWO Positioning Aids



...realising optimum x-ray image performance



CAWO Positioning Aids

Help the radiographer in setting-up the patient for any particular examination and therapy. They are made from radiotransparent polyester foam. The skin friendly polyurethane cover is dirt-resistant and easy to clean.

Available in a variety of shapes, suitable for different set-ups.

No. 1:	Rectangular Pad 250 x 200 x 70 mm	No. 7:	Roll 250 x 100 mm
No. 2:	Large Wedge 250 x 200 x 110 mm	No. 8:	A-P Head Mould 250 x 200 x 100 mm
No. 3:	Small Wedge 120 x 200 x 60 mm	No. 9:	P-A- Head Mould 250 x 200 x 100 mm
No. 4:	Flat wedge 250 x 200 x 700 mm	No. 10:	Right Lateral Head Mould 250 x 250 x 100 mm
No. 5:	Round Pad 200 x 50 mm	No. 11:	Left Lateral Head Mould 250 x 250 x 100 mm
No. 6:	Ring Pad 200 x 50 mm	No. 12:	Baby Mould 400 x 300 x 120 mm
		No. 13:	Child Mould 340 x 540 x 120mm

Printed in November 2006

Pads with robust and non-irritant Nylon:

2000 x 563 x 10 mm,
2000 x 563 x 20 mm,
2000 x 563 x 30 mm

Standard set:

Standard set with box:

2 x No. 2, No. 3 and No. 4 / 1 x No. 1, No. 5, No. 6 and No. 7
2 x No. 2, No. 3 and No. 4 / 1 x No. 1, No. 5, No. 6 and No. 7

Cleaning: with water and currently available cleaning agent. avoid contact with alcohol or trichlorine-ethylene.

CAWO Photochemisches Werk GmbH
P.O. Box 1129


D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-0

Fax: (++49) 8252 9109-19

Internet: www.cawo.com

E-mail: cawo@cawo.com



Characteristics of the products described in this publication can be changed at any time without notice

