

MR-CONDITIONAL STRETCHER



This heavy-duty stretcher is an ideal unit for MRI patient transport. Finished in attractive white enamel, it comes complete with mattress, restraining straps, and retractable guard rails. Although light in weight for easy, controlled mobility, the MRI Stretcher holds patients safely and comfortably. It rolls on 5" ball-bearing swivel casters that lock firmly to secure the stretcher in place when positioned.

MR-Conditional Stretcher

- Load capacity: 500lbs (227 kg)
- Weight: 85lbs (30 kg)
- Dimensions: 72"L x 25"W x 33"H (183 x 64 x 83 cm)

Model: MRST, USA