

## Non-Magnetic Laryngoscope

### Features:

- Non-Magnetic Fiber Optic Handle
- Low Magnetic Field Batteries
- Full Assortment of Blades

### Available Non-Magnetic Blades:

- AM MAC: 0, 1, 2, 3, 3.5, 4
- ENG MAC: 1, 2, 3, 4
- MILLER: 0, 1, 2, 3, 4
- Wisconsin: 2, 3, 4

### Safety Note:

- Do Not use these blades in any handle that is not marked "MRI TESTED"
- The Fiber Optic Laryngoscope Handle was tested safe to 1.5 Tesla static magnetic field.



### Item Number:

### Description:

|           |                                                              |
|-----------|--------------------------------------------------------------|
| 03-310001 | Non-Magnetic Laryngoscope handle - battery and bulb included |
| 03-310010 | Low Magnetic Field Batteries 2/pk                            |
| 03-310011 | Pk Bulbs 2/pk                                                |