



TRLYA - MR Conditional Adjustable Height Trolley

Tested for use at 7 Tesla, these Trendelenberg trolleys feature an adjustable back rest – gas spring assisted, IV pole with 4 fixing points, 4 low rolling resistance casters, all with foot operated brakes. Corner Buffering for protection. Folding Cot Sides. Flame retardant, high density mattress in a wipe clean material to prevent the spread of infection.

- Weight Capacity 220 Kg (490lbs)
- Base Size 1475mm (58") x 610mm (23½") Maximum Height 940mm (37")
- Minimum Height 495mm (19½") To table top.
- Mattress Thickness 50mm (2")
- Overall Top size 2175mm (85½") x 640mm (25¼")