

MRI Stretcher

This heavy-duty, all aluminum non-magnetic stretcher is an ideal unit for MRI patient transport. The stretcher has a welded frame and is finished in white powder coat. It comes complete with mattress, restraining straps and retractable guard rails. Although light in weight for easy, controlled mobility, the MRI Stretcher accommodates up to a 500-lb patient safely and comfortably. It rolls on 5" casters for easy maneuverability.

Specifications:

03-PH1004 and 03-PH1005

- Dimensions:
 - Tabletop: 72"L x 24"W (182.8 x 60.9 cm)
 - Overall: 74.5"L x 27.75"W x 34"H (189.3 x 70 x 86 cm)
- IV Pole Socket (iv pole ordered separately)
- Wheels: 5" locking casters (independently locking)
- Finish: Powder coat
- Patient Capacity: 500 lb (226.8 kg);
weight tested to four times the patient load rating.
- Compliance: MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System.
- Weight: 85 lb (39kg) (03-PH1004 only)



03-PH1005 only

- Fowler Back: 0° to 80° infinitely adjustable, manual
- Weight: 110 lb (49.9kg)

Item Number:

Description:

03-PH1004

MRI Fixed Stretcher

03-PH1005

MRI Fixed Stretcher with Elevating Headrest (Fowler Back)

03-PHPOLE

IV Pole for Stretcher

CT-05-11, CORPORATE TOWER SUBANG SQUARE, JALAN SS15/4G, SUBANG JAYA, 47500, SELANGOR DARUL ESHAN, MALAYSIA

TEL:603-56366732

FAX:603-56367732

Web Page: www.hnamedical.com