


MRI Patient Deluxe Positioning Kit

These positioning products offer cost effective solutions for ensuring correct and comfortable positioning during MRI procedures. Proper patient positioning is critical to accurate scans, diagnosis, and helps to minimize the cost of each scan.



Standard Foam

- Open Cell Foam
- Radiolucent
- Fire Retardent
- Economic General Positioning

Vinyl Cover (-20)

- Extremely Durable Vinyl Cover for standard fom positioning devices
- Heat Sealed to prevent fluid absorption/contamination

(not recommended for low KVP techniques)

Scan Coat (-SCB)

- Strong, State-of-the-Art Sprayed Vinyl Coating:
Withstands normal clinical use without tearing or peeling
- Very Low Attenuation
- Skid Resistant
- Resists Contamination

Item Number: **Description:**

03-PK-1-20 Deluxe Positioning Kit - Vinyl Coated

Description:	Quantity
Wedges 3"H x 3"W x 7.25"L	2
Wedge 6.75"H x 6.75"W x 7.25"L	1
Circular disc 1.5"H x 7"D	1
Wedges 2.25"H x 9.5"W x 7.25"L	2
Rectangle 2"H x 7.5"W x 9.5"L	1
Decubitus block 4"H x 18"W x 24"L	1
Cylinder 12"H x 4.25"D	1
Thin mattresses .25"H x 18"W x 24"L	2
Table pad 1"H x 24"W x 72"L	1