

MRI MOBILE TRANSPORTER COT



MRI Mobile Transporter Cot

The Mobile Transporter, Multi-Level Roll-In Cot combines the benefits of a multi-level cot with the easy loading of a roll-in style cot. Empty one-person load capability for user convenience. Entirely non-magnetic for use in MRI facilities. Utilizes a manual backrest. Includes burgundy bolster mattress and three black patient restraints.

- Lightest cot available with a 500-pound load capacity
- Sturdy X-frame undercarriage is easily raised by second EMT during loading
- Standard with 4 locking casters
- Dual ratchet bar adjustment mechanism is simple to operate and permits eight levels of height adjustment
- Safety U-bar and hook holds the Mobile Transporter at the rear of the ambulance during unloading as undercarriage is lowered to the ground
- Low foot-end lifting frame enables shorter EMTs to get better leverage when lifting the Mobile Transporter
- Head-end drop frame and foot-end lead handle for improved control/maneuverability
- Adjustable backrest angle from 0° to 75° allows patient to be in a comfortable position
- Built-in shock position permits instant change in patient positioning by elevating feet 15 degrees (approximately 10") without compromising patient's breathing
- Improved lubrication-free easy-roll wheels for greater cot mobility
- MR-conditional to 3-Tesla

Newmatic , MCMTC