

MRI Non-Magnetic 20" Wide Wheelchair w/ IV Pole and Cylinder Holder



This 20" MRI non-magnetic wheelchair is constructed entirely of non-ferrous materials. Specifically designed for use in and around the MRI suite. The storage pocket on rear of chair has "MR" prominently displayed, reducing the chance it will be confused with a standard magnetic wheelchair or end up in another department.

Specifications:

- Seat Width: 20"
- Overall Width: 28"
- Overall Depth: 44" w/ footrests
- Seat Depth: 16"
- Overall Height: 38"
- Back Height: 18"
- Front Casters: 8" x 1" solid
- Rear Wheel: 24" solid
- Wheel Locks: push to lock
- Arm Rests: Removable
- Foot Rests: Removable
- Upholstery: Burgundy
- Weight: 46 lbs.
- MR Conditional to: 3-Tesla
- Weight Capacity: 250 lbs.