

MR-CONDITIONAL LED LARYNGOSCOPE KIT



MR-Conditional LED Laryngoscope Kit

- Tested and demonstrated to be MR-conditional for magnetic field interactions at 3-Tesla with highest spatial gradient field of 720 gauss/cm or less
- LED light is cool running long lasting and uses very little power to operate
- Non-magnetic steel autoclavable fiberoptic blades and light guide
- Powered by low magnetic battery for safety in the MRI suite
- Distinctly identifiable silver blades with gold tubes and black fiber optic handles reduce risk of incompatible instruments getting into MRI suites
- Laser printed blades and handles available in a pre-packaged set or sold individually
- Special material blades are developed and approved for MR environments
- Blades available in a variety of sizes in English Mac, American Mac, Miller, and Wisconsin styles

MRI-Conditional LED Laryngoscope Kit

- 6 MR-conditional blades (Macintosh blades 2, 3, 4 & Miller blades 1, 2, 3)
- 2 MR-conditional special battery packs
- 1 bulb, pre-installed in handle
- 1 MR-conditional medium-size handle
- Nonferrous storage case