

MRI COMFORT STOOL WITH GAS LIFT



Model: CSTL



Model: CSTLB



Model: CSTLBA

The most comfortable stool you'll ever sit on even for extended periods of time. This ultra-cushioned MRI stool features a non-magnetic gas lift the first of its kind. It sits on an ultra-stable 5-point base with easy-glide swivel casters for easy mobility. Black vinyl upholstery. Tested MR-conditional to 7-Tesla.

MRI Comfort Stool with Gas Lift

- Smooth-rise non-magnetic gas lift
- 4" thick seat cushion; 16"Dia
- Seat height range from 19"-24"
- MR-conditional to 7-Tesla
- Five-star base with easy-glide swivel casters

Model:CSTL: MRI Comfort Stool with Gas Lift

Model: CSTLB: MRI Comfort Stool with Gas Lift and Back

Model: CSTKBA: MRI Comfort Stool with Gas Lift, Back, and Arms