



Coarto

DIGITAL MAMMOGRAPHY FORCE GAUGE

Coarto is a digital force gauge which allows the user to measure breast compression on mammography systems.

Coarto meets the following standards:

IPEM Report 89
EUREF 2a.2.1.4

Coarto includes the following:

Breast compression scale
-maximum measurement 300N (or 30kg)
-accurate to 5g
-dynamic display
-peak hold display

Compressible rubber block



For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102