



# PIXMAM

## DIGITAL MAMMOGRAPHY

PIXMAM is a digital mammography image quality phantom comprising a stack of finely tolerated PMMA plates which meets the routine testing requirements of the following standards;

- NHSBSP pub 604-3 :2009
- NHSBSP pub 702-1 :2007
- EUREF 4th edition of the European Guidelines for Quality Assurance in Breast Cancer Screening and Diagnosis

PIXMAM comprises PMMA plates of dimensions;

(1pc) 300 x 240 x 20mm  
(including encapsulated 99.99% Al foil)  
(2pc) 300 x 240 x 20mm  
(1pc) 300 x 240 x 10mm  
(1pc) 300 x 240 x 5mm

PIXMAM allows the following tests to be conducted;

- AEC checks
- detector homogeneity
- SNR
- CNR



PIXMAM 300 x 240



PIXMAM 250d

For more information or to find your regional distributor please visit our website

[www.leadstestobjects.com](http://www.leadstestobjects.com) | tel. +44(0)1423 321102



# PIXMAM

DIGITAL MAMMOGRAPHY



PIXMAM RADIOGRAPH

For more information or to find your regional distributor please visit our website

[www.leedstestobjects.com](http://www.leedstestobjects.com) | tel. +44(0)1423 321102

Copyright © 2011 Leeds Test Objects Ltd | SpecSh 275c 16/02/2012