

SEDENTEXCT IQ MINI

DENTAL CBCT

PMMA cylinder with;

- 160 mm diameter
- 144 mm length
- density $1.20 \pm 0.01 \text{ g/cm}^3$

The body of the cylinder allows the following test to be undertaken;

Noise/Uniformity

- The upper and lower sections of the phantom are uniform PMMA (density $1.20 \pm 1.00\%$)

The test feature in the centre of the cylinder is included to perform the following measurement;

Pixel Intensity Value

- 25 mm diameter discs of Aluminium, PTFE, delrin and LDPE suspended in PMMA



SEDENTEXCT IQ MINI

For more information or to find your regional distributor please visit our website