

SEDENTEXCT SWIFT

DENTAL CBCT

Automatic scoring software for the SedentexCT IQ phantom

SedentexCT SWIFT includes the following tools;

Selecting tool

- enabling selection and evaluation of the inserts



Zooming tool

- enabling zooming in and out of the image on the screen



Measuring tool

- enabling distance measurements between two points (in mm)



Window level tool

- enabling adjustments in brightness / contrast




fig. 1, selecting tool




fig. 2, measuring tool

For more information or to find your regional distributor please visit our website