



TO BIOPSY

MAMMOGRAPHY

Leeds Test Object TO BIOPSY is designed to be used to check the localisation performance of stereotactic core needle biopsy systems used in mammography.

TO BIOPSY consists of a PMMA block containing two sets of holes at different depths and diameters.

TO BIOPSY enables the following tests to be made:

- Fine needle aspiration stereotactic localisation testing
- Core biopsy stereotactic localisation testing



For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102