



TO DR

DYNAMIC RANGE

These test objects are designed for systems which operate by sensing mean signal levels, the TO DR tests dynamic range by measuring detail visibility over a wide range of input Xray intensities. Different designs cater for subtraction and Non-subtraction systems.

The TO DR (4 piece)

Supplied in a set that allows it to be used for both subtraction and non-subtraction systems. Comprises of the following test objects:

- Filter Plate
- Non-Subtraction Detail Plate
- Subtraction Detail Plate
- Blank Plate



The TO DR (2 piece)

For use with non-subtraction systems only. Comprises of the following test objects:

- Filter Plate
- Non-Subtraction Detail Plate

For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102