



TO FSJ Mammo

MAMMOGRAPHY

The focal spot measurement jig consists of a cassette holder which can be set at any angle up to 20° above the horizontal.

Mounted on the cassette holder is a tower which holds a 1° star pattern, supplied with the jig, giving a magnification factor approximately 2.5 x for an F.F.D. of 600mm.

At the base of the tower is a centre-marker (a circular ring) which, when concentric with the star pattern image, indicates that the jig is correctly aligned to the central beam axis.



For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102