

X-ray viewers without shutters

**For film selection, for checking the
adjustment technique, for demonstration**



MOD

**the favorably priced module system with
individually controlled viewing sections**

BW I, II

**viewing panels in special designs made to
customer specifications, with 1 or 2 viewing
areas, one arranged on top of the other**

The favorably priced Modul system with individually controlled viewing sections for film selection and for the checking of the adjustment technique

- ¹⁾ With reference to the fluorescent lamp luminous density indicated in cd/m², allowance has been made for a light output loss of roughly 10 per cent.

Description								Order No
Smooth electronic brightness control, adjustment range approx. 50%, each viewing section can be individually controlled								201005
1 pair of viewbox legs (for desk model)								301024

Type	Viewing area in cm W x H	Housing dimensions in mm W H D			Number and wattage of fluoresc. lamps	Power input in watt	Current in ampere	Weight in kgs (approx.)	Order No
MOD	35 x 43	374	520	135	2 x 15 W	40	0,4	9	103040
	70 x 43	724	520	135	4 x 15 W	80	0,8	18	103041
	105 x 43	1074	520	135	6 x 15 W	120	1,2	27	103042
	140 x 43	1424	520	135	8 x 15 W	160	1,6	36	103043

